

| VALIDATION RESULTS FOR METHOD NO. Covance Harrogate 2001-032-D1 | | | | | | |
|---|-------------------------------|-----------------|---|-----------|------------|--|
| ANALYTE | MDA | MEAN | RELATIVE STD. DEV.(%) | RANGE | SEE (conc) | |
| MATRIX | Human plasma | | | | | |
| DATE | 15-Aug-01 | | | | | |
| LLOQ | (Provide detection criterion) | 0.4 µmol/L | | | | |
| | | Accuracy 107.6% | | | | |
| | | RSD 6.0% | | | | |
| STANDARD ERROR OF ESTIMATE OF QC SAMPLES (concentration units) | | | | | | |
| (SEE = S.D./SQRT (# of determinations)) | | | | | | |
| | 1 µmol/L | Low QC | | | 0.012 | |
| | 4 µmol/L | Med QC | | | 0.034 | |
| | 8 µmol/L | High QC | | | 0.068 | |
| DESCRIPTIVE INFORMATION ON INSTRUMENT RESPONSE | | | | | | |
| Function (Linear, Quadratic, etc.) | Linear | | | | | |
| Slope (mean) | 67927 | | | | | |
| Intercept (mean) | 7935 | | | | | |
| Other (mean regression value) | 0.98855(1/x weighted) | | | | | |
| RECOVERY (as a %) | | | | | | |
| High QC in (matrix)= | 8 µmol/L | 92.5% | Calculated as a % of the peak area observed for derivatised water standards | | | |
| Med QC in (matrix)= | 4 µmol/L | 94.7% | | | | |
| Low QC in (matrix)= | 1 µmol/L | 100.0% | | | | |
| Results with different sample amounts (results of QC evaluations) intra-assay | | | | | | |
| LLOQ QC= | 0.4 µmol/L | 0.43 | 6.0% | 0.39-0.45 | 0.013 | |
| Low QC= | 1 µmol/L | 0.94 | 4.1% | 0.90-1.00 | 0.016 | |
| Med QC= | 4 µmol/L | 3.90 | 2.2% | 3.79-4.05 | 0.035 | |
| High QC= | 8 µmol/L | 7.55 | 2.4% | 7.34-7.82 | 0.075 | |
| Results with different sample amounts (results of QC evaluations) inter-assay | | | | | | |
| LLOQ QC= | n/a | n/a | n/a | n/a | n/a | |
| Low QC= | 1 µmol/L | 0.95 | 6.1% | 0.84-1.10 | 0.012 | |
| Med QC= | 4 µmol/L | 3.94 | 4.3% | 3.54-4.15 | 0.034 | |
| High QC= | 8 µmol/L | 7.50 | 4.5% | 6.85-8.10 | 0.068 | |
| STABILITY (assuming any degradation is linear): | | | | | | |
| a) Room Temperature in Matrix (+/- %/12 hrs) | | 7.7% maximum | | | | |
| b) Freeze/thaw (+/- %/3 cycles) | | 3.3% maximum | | | | |
| c) Processed sample (+/- %/24 hrs) | TBA | % maximum | | | | |

| Headings | | | | |
|------------------|--------|--------|--------|--|
| 1A | M10A | Std 10 | 44.000 | |
| Units | mg/dl | Cal 1 | 0.4 | |
| Species | human | Cal 2 | 0.6 | |
| Matrix | plasma | Cal 3 | 1 | |
| Height area | area | Cal 4 | 2 | |
| DP% | 1 | Cal 5 | | |
| DP/area compound | 1 | Cal 6 | 3 | |
| ULO QC conc | 0.4 | Cal 7 | 7.5 | |
| Lo QC conc | 1 | Cal 8 | 9 | |
| Mid QC conc | 4 | Cal 9 | 10 | |
| Hi QC conc | 8 | | | |

Footnotes from SOP 66 Rev 1: Batch annotation

AD = Above data system

AR = Anomalous result

BD = Below data system

CO = Confirms original result

DN = Data not used

DNR = Data not required

HI = High internal standard

M = Instrument malfunction (add reason if known)

P = Interferent peak

S = Insufficient sample

L = Low internal standard

A = This result is a median value

ME = Mis-extraction

MI = Mis-injection

NR = No result

NS = No sample taken/received

OR = Over range

PC = Poor chromatography

PF = Poor fit with rest of regression (excluding standard improves the correlation coefficient / curve fit)

PT = Poor fit, calibration standard outside acceptance criteria

PP = Positive predose

PS = Poor sensitivity

SA = Sample dropped/spilt/spoilt during preparation (Lab accident)

SE = Spiking error

UD = Under range (when LLOQ not met)

UN = Unknown result (when not a standard or a sample)